

Dry Active series

HIGH TECHNOLOGY ELECTRONIC
DRY CABINET



DRI Active

FEATURES

Stores various products sensitive to humidity change: Camera, Electronic goods, Reagent, Specimen, Medicine, Semi-conductor etc.

High-density foam gaskets perfectly keep contact from the outside

It is lighter than that of other companies, and excellent for protecting fragileness as all windows are made of PC

A wide digital LCD display screen indicates internal humidity and temperature.

Auto-Dry Cabinet devices stably and fastly keep internal humidity at 20%-40% (25%-40% in case of Aliminum Desiccator Cabinet) by adopting electronic dehumidifying methods.

Dry Active-Desiccator Cabinet



Catalog no:	Description	Sales Unit
KA.33-70	General Type	EA
KA.33-80	Shelf (option)	EA

■ Material :

Body : ABS + PC

Frame : ABS (Acrylonitrile Butadiene Styrene)

■ Size : 47.6L

Outside Dimensions : 300(W)×380(D)×520(H)mm

Inside Dimensions : 280(W)×340(D)×500(H)mm

Size of shelf : 280(W)×340(D)mm

■ Accessories :

3 Shelves, 1 Desiccant Tray and 1 Digital Thermo-Hygrometer(built-in) included

■ Digital Thermo & Hygrometer Specifications :

Indicating humidity : 10-99%

Deviation of humidity : ±5%

Indicating temperature : -50°C~70°C (-58°F~158°F)

Deviation of temperature : ±1°C(±2°F)

Auto Desiccator Cabine (Dry Active)



Catalog no:	Description	Sales Unit
KA.33-71	Auto Type	EA
KA.33-80	Shelf (option)	EA

Material :

Body : ABS + PC

Frame : ABS(Acrylonitrile Butadiene Styrene)

■ Size : 47.6L

Outside Dimensions: 300(W)×380(D)×635(H)mm Inside Dimensions: 280(W)×340(D)×500(H)mm Size of shelf: 280(W)×340(D)mm

■ Accessories :

3 Shelves, 1 Digital Thermo-Hygrometer(built-in) included

■ Digital Thermo & Hygrometer Specifications : Indicating humidity : 10-99%

Deviation of humidity : ±5%

Indicating temperature : -50°C~70°C (-58°F~158°F) Deviation of temperature : ±1°C(±2°F)

Auto Desiccator Cabine (Dry Active)



Catalog no:	Description	Sales Unit
KA.33-73	Auto-Humidity Controller Type	EA
KA.33-80	Shelf (option)	EA

■ Material :

Body : ABS + PC

Frame : ABS(Acrylonitrile Butadiene Styrene)

■ Size : 47.6L

Outside Dimensions : 300(W)×380(D)×635(H)mm Inside Dimensions : 280(W)×340(D)×500(H)mm Size of shelf : 280(W)×340(D)mm

■ Accessories :

3 Shelves, 1 Digital Thermo-Hygrometer(built-in) included

■ Humidity Control Range : 20-40%

■ Digital Thermo & Hygrometer Specifications :

Indicating humidity : 10-80%

Deviation of humidity : ±3%

Indicating temperature : -50°C~99°C (-58°F~210°F) Deviation of temperature : ±1°C (±2°F)

■ Power Supply : AC 220V, 50/60Hz

Desiccator Cabinet (Dry Active) – UV Protection



Catalog no:	Description	Sales Unit
KA.33-70X	General Type - Amber Color	EA
KA.33-80	Shelf (option)	EA

Material :

Body : ABS + PMMA

Frame : ABS(Acrylonitrile Butadiene Styrene)

■ Size : 47.6L

Outside Dimensions: 300(W)×380(D)×520(H)mm

Inside Dimensions: 280(W)×340(D)×500(H)mm

Size of shelf: 280(W)×340(D)mm

■ Accessories :

3 Shelves, 1 Desiccant Tray and 1 Digital Thermo-Hygrometer(built-in) included

■ Digital Thermo & Hygrometer Specifications :

Indicating humidity : 10-99%

Deviation of humidity : ±5%

Indicating temperature : -50°C~70°C (-58°F~158°F)

Auto Desiccator Cabine (Dry Active)



Catalog no	Description	Sales Unit
KA.33-74	Auto-Humidity Controller Type	EA
KA.33-80	Shelf (option)	EA

■ Material :

Body : ABS + PC

Frame : ABS(Acrylonitrile Butadiene Styrene)

■ Size : 95.2L

Outside Dimensions : 300(W)×380(D)×1130(H)mm Inside Dimensions : 280(W)×340(D)×1000(H)mm Size of shelf : 280(W)×340(D)mm

■ Accessories :

6 Shelves, 1 Digital Thermo-Hygrometer(built-in) included

■ Humidity Control Range : 20-40%

■ Digital Thermo & Hygrometer Specifications :

Indicating humidity : 10-80%

Deviation of humidity : ±3%

Indicating temperature : -50°C~99°C(-58°F~210°F) Deviation of temperature : ±1°C(±2°F)

■ Power Supply : AC 220V, 50/60Hz

Auto Desiccator Cabinet (Dry Active) – UV Protection



Catalog no:	Description	Sales Unit
KA.33-73X	Auto-Humidity Controller Type - Amber Color	EA
KA.33-80	Shelf (option)	EA

Material :

Body : ABS + PMMA

Frame : ABS(Acrylonitrile Butadiene Styrene)

■ Size : 47.6L

Outside Dimensions: 300(W)×380(D)×635(H)mm

Inside Dimensions: 280(W)×340(D)×500(H)mm

Size of shelf: 280(W)×340(D)mm

■ Accessories :

3 Shelves, 1 Digital Thermo-Hygrometer(built-in) included

■ Humidity Control Range : 20-40%

■ Digital Thermo & Hygrometer Specifications :

Indicating humidity : 1-99%

Deviation of humidity : ±3%

Indicating temperature : -50°C~99°C (-58°F~210°F)

Deviation of temperature : ±1°C(±2°F)

■ Power Supply : AC220V, 50/60Hz

Auto Desiccator Cabinet (Dry Active) – UV Protection



Catalog no:	Description	Sales Unit
KA.33-74X	Auto-Humidity Controller Type - Amber Color	EA
KA.33-80	Shelf (option)	EA

Material :

Body : ABS + PMMA

Frame : ABS(Acrylonitrile Butadiene Styrene)

■ Size : 95.2L

Outside Dimensions: 300(W)×380(D)×1130(H)mm

Inside Dimensions: 280(W)×340(D)×1000(H)mm

Size of shelf: 280(W)×340(D)mm

■ Accessories : 6 Shelves, 1 Digital Thermo-Hygrometer(built-in) included

■ Humidity Control Range : 25~40%

■ Digital Thermo & Hygrometer Specifications :

Indicating humidity : 1-99%

Deviation of humidity : ±3%

Indicating temperature : -50°C~99°C (-58°F~210°F)

Deviation of temperature : ±1°C(±2°F)

■ Power Supply : AC220V, 50/60Hz

Aluminum Desiccator Cabinet – UV Protection



Catalog no:	Description	Sales Unit
KA.33-76X	500x525x1250mm(General Type) - Amber Color	EA
KA.33-76AX	500x525x1250mm(Auto Type) - Amber Color	EA
KA.33-84	Shelf (Option)	EA

■ Material :

Body : PC / Frame : Aluminum / Shelf : Stainless steel

■ Size : 203.70L

Inside Dimensions : 420 (W)×485 (D)×1000(H)mm

Size of shelf : 416(W)×476(D)mm

■ Accessories : 3 Shelves, 1 Desiccant Tray(Only for General Type), 1 Digital Thermo-Hygrometer(built-in)

■ Humidity Control Range : 25~40%

■ Digital Hygrometer Specifications :

Indicating humidity : 10~99%

Deviation of humidity : ±5%

Indicating temperature : -50°C~70°C(-58°F~158°F)

Deviation of temperature : ±1°C (±2°F)

■ Power Supply : AC 220V, 50/60Hz

Aluminum Desiccator Cabinet – UV Protection



Catalog no:	Description	Sales Unit
KA.33-75X	500x525x850mm(General Type) - Amber Color	EA
KA.33-75AX	500x525x1060mm (Auto Type) - Amber Color	EA
KA.33-84	Shelf (option)	EA

■ Material :

Body : PC / Frame : Aliminum / Shelf : Stainless steel

■ Size : 167.03L

Inside Dimensions : 420 (W)×485 (D)×820(H)mm

Size of shelf : 416(W)×476(D)mm

■ Accessories :2 Shelves, 1 Desiccant Tray(Only for General Type), Digital Thermo-Hygrometer(built-in)

■ Humidity Control Range : 25-40%

■ Digital Hygrometer Specifications :

Indicating humidity : 10-99%

Deviation of humidity : ±5%

Indicating temperature : -50°C~70°C(-58°F~158°F)

Deviation of temperature : ±1°C (±2°F)

■ Power Supply : AC 220V, 50/60Hz

Aliminum Desiccator Cabinet – UV Protection



Catalog no:	Description	Sales Unit
KA.33-77X	500x525x1900mm(General Type) - Amber Color	EA
KA.33-77AX	500x525x1900mm(Auto Type) - Amber Color	EA
KA.33-84	Self(option)	EA

■ Material :

Body : PC / Frame : Aliminum / Shelf : Stainless steel

■ Size : 338.14L

Inside Dimensions : 420 (W)×485 (D)×1660(H)mm

Size of shelf : 416(W)×476(D)mm

■ Accessories :

4 Shelves, 1 Desiccant Tray(Only for General Type), 1 Digital Thermo-Hygrometer(built-in)

■ Humidity Control Range : 25-40%

■ Digital Hygrometer Specifications :

Indicating humidity : 0-100%

Deviation of humidity : ±3%

■ Power Supply : AC 220V, 50/60Hz

Auto Desiccator Cabinet (Dry Active) – UV Protection



Catalog no:	Description	Sales Unit
KA.33-74X	1010x525x1900mm(General Type)	EA
KA.33-80	1010x525x1900mm(Auto Type)	EA
KA.33-84	Shelf (option)	EA

■ Material :

Body : PC / Frame : Aliminum / Shelf : Stainless steel

■ Size : 748.74L

Inside Dimensions : 930 (W)×485 (D)×1660(H)mm

Size of shelf : 416(W)×476(D)mm

■ Accessories :

8 Shelves, 2 Desiccant Tray(Only for General Type), Digital Hygrometer(built-in)

■ Humidity Control Range : 25-40%

■ Digital Hygrometer Specifications :

Indicating humidity : 0-100%

Deviation of humidity : ±3%

■ Power Supply : AC 220V, 50/60Hz

Aliminum Desiccator Cabinet



Catalog no:	Description	Sales Unit
KA.33-76	500x525x1250mm(General Type)	EA
KA.33-76A	500x525x1250mm(AutoType)	EA
KA.33-84	Shelf (Option)	EA

■ Material :

Body : PC / Frame : Aliminum / Shelf : Stainless steel

■ Size : 203.70L

Inside Dimensions : 420 (W)×485 (D)×1000(H)mm

Size of shelf : 416(W)×476(D)mm

■ Accessories :3 Shelves, 1 Desiccant Tray(Only for General Type), 1 Digital Hygrometer(built-in)

■ Humidity Control Range : 25-40%

■ Digital Hygrometer Specifications :

Indicating humidity : 10-99%

Deviation of humidity : ±5%

Indicating temperature : -50°C~70°C(-58°F~158°F)

Deviation of temperature : ±1°C (±2°F)

■ Power Supply : AC 220V, 50/60Hz

Aluminum Desiccator Cabinet



Catalog no:	Description	Sales Unit
KA.33-75	500x525x850mm(General Type) - Amber Color	EA
KA.33-75AX	500x525x1060mm (Auto Type) - Amber Color	EA
KA.33-84	Shelf (option)	EA

■ Material :

Body : PC / Frame : Aliminum / Shelf : Stainless steel

■ Size : 167.03L

Inside Dimensions : 420 (W)×485 (D)×820(H)mm

Size of shelf : 416(W)×476(D)mm

■ Accessories : 2 Shelves, 1 Desiccant Tray(Only for General Type), Digital Hygrometer(built-in)

■ Humidity Control Range : 25-40%

■ Digital Hygrometer Specifications :

Indicating humidity : 10-99%

Deviation of humidity : ±5%

Indicating temperature : -50°C~70°C(-58°F~158°F)

Deviation of temperature : ±1°C (±2°F)

■ Power Supply : AC 220V, 50/60Hz

Aliminum Desiccator Cabinet



Catalog no:	Description	Sales Unit
KA.33-77	500x525x1900mm(General Type)	EA
KA.33-77A	500x525x1900mm(Auto Type)	EA
KA.33-84	Self(option)	EA

■ Material :

Body : PC / Frame : Aliminum / Shelf : Stainless steel

■ Size : 338.14L

Inside Dimensions : 420 (W)×485 (D)×1660(H)mm

Size of shelf : 416(W)×476(D)mm

■ Accessories :

4 Shelves, 1 Desiccant Tray(Only for General Type), 1 Digital Hygrometer(built-in)

■ Humidity Control Range : 25-40%

■ Digital Hygrometer Specifications :

Indicating humidity : 0-100%

Deviation of humidity : ±3%

■ Power Supply : AC 220V, 50/60Hz

Auto Desiccator Cabinet



Catalog no:	Description	Sales Unit
KA.33-78	1010x525x1900mm(General Type)	EA
KA.33-78A	1010x525x1900mm(Auto Type)	EA
KA.33-85	Shelf (option)	EA

■ Material :

Body : PC / Frame : Aliminum / Shelf : Stainless steel

■ Size : 748.74L

Inside Dimensions : 930 (W)×485 (D)×1660(H)mm

Size of shelf : 430(W)×476(D)mm

■ Accessories :

8 Shelves, 2 Desiccant Tray(Only for General Type), Digital Hygrometer(built-in)

■ Humidity Control Range : 25-40%

■ Digital Hygrometer Specifications :

Indicating humidity : 0-100%

Deviation of humidity : ±3%

■ Power Supply : AC 220V, 50/60Hz

Gas Exchangeable Desiccator Cabinet



Catalog no:	Description	Sales Unit
KA.33-76GE	500x525x1250mm (Gas Exchangeable type)	EA

■ Material :

Body : PC / Frame : Aliminum / Shelf : Stainless steel

■ Size : 203.70L

Inside Dimensions : 420 (W)×485 (D)×1000(H)mm

Size of shelf : 416(W)×476(D)mm

■ Accessories :

3 Shelves, 1 Desiccant Tray(Only for General Type), Digital Hygrometer(built-in)

■ Digital Hygrometer Specifications :

Indicating humidity : 10-99%

Deviation of humidity : ±5%

Indicating temperature : -50°C~70°C(-58°F~158°F)

Deviation of temperature : ±1°C (±2°F)

■ Gas in/out Valve Diameter Ø 6mm (Ø 6mm Tube Recommended)

Gas Exchangeable Desiccator Cabinet

Catalog no:	Description	Sales Unit
KA.33-75GE	500x525x850mm(Gas Exchangeable type)	EA



- Material :
Body : PC / Frame : Aliminum / Shelf : Stainless steel
- Size : 167.03L
Inside Dimensions : 420 (W)×485 (D)×820(H)mm
Size of shelf : 416(W)×476(D)mm
- Accessories :
2 Shelves, 1 Desiccant Tray(Only for General Type), Digital Hygrometer(built-in)
- Digital Hygrometer Specifications :
Indicating humidity : 10-99%
Deviation of humidity : ±5%
Indicating temperature : -50°C~70°C(-58°F~158°F)
Deviation of temperature : ±1°C (±2°F)
- Gas in/out Valve Diameter Ø 6mm (Ø 6mm Tube Recommended)

Gas Exchangeable Desiccator Cabinet

Catalog no:	Description	Sales Unit
KA.33-77GE	500x525x1900mm(Gas Exchangeable type)	EA



- Material :
Body : PC / Frame : Aliminum / Shelf : Stainless steel
- Size : 338.14L
Inside Dimensions : 420 (W)×485 (D)×1660(H)mm
Size of shelf : 416(W)×476(D)mm
- Accessories :
4 Shelves, 1 Desiccant Tray(Only for General Type), Digital Hygrometer(built-in)
- Digital Hygrometer Specifications :
Indicating humidity : 10-99%
Deviation of humidity : ±5%
Indicating temperature : -50°C~70°C(-58°F~158°F)
Deviation of temperature : ±1°C (±2°F)
- Gas in/out Valve Diameter Ø 6mm (Ø 6mm Tube Recommended)